


A35 NITRILE EXAMINATION GLOVES



FEATURES	SPECIFICATIONS		PACKAGING		
<ul style="list-style-type: none">• Powder-Free• Ambidextrous• Beaded Cuff• Non-Sterile• Not Made From Natural Rubber Latex	MATERIAL:	Nitrile	XS-L: 100 gloves/box 10 boxes/case		
	POWDER-FREE:	Yes			
	COLOR:	Blue	XL: 90 gloves/box 10 boxes/case		
	GLOVE LENGTH:	245mm / 9.5 inches			
	FINGERTIP SURFACE:	Textured			
	PALM THICKNESS:	0.08mm / 3.1mil			
	FINGER THICKNESS:	0.11mm / 4.3mil			
	AVAILABLE SIZES:	XS, S, M, L, XL			
	INSPECTION LEVEL (FFH):	AQL 1.5			
	NOT MADE FROM NATURAL RUBBER LATEX:	Yes			
RECOMMENDED APPLICATIONS			ORDERING INFORMATION		
<ul style="list-style-type: none">• Medical Examinations• Laboratory Research• Food Handling• Janitorial Work• Light Duty Maintenance			ITEM	SIZE	QUANTITY
			A350	XS	100/Box, 1000/Case
			A351	S	100/Box, 1000/Case
			A352	M	100/Box, 1000/Case
			A353	L	100/Box, 1000/Case
			A354	XL	90/Box, 900/Case

CHEMICAL PERMEATION TESTING PER ASTM D6978:

RESULTS SHOWN FOR EACH AND EVERY CHEMICAL LISTED:
AVERAGE BREAKTHROUGH DETECTION TIME: >240 MINUTES
ADDITIONAL DATA AVAILABLE UPON REQUEST

CHEMICALS TESTED

CISPLATIN 1.0 MG/ML (1,000PPM)	ETOPOSIDE 20.0 MG/ML (20,000PPM)	PACLITAXEL 6.0 MG/ML (6,000PPM)
CYCLOPHOSPHAMIDE 20.0 MG/ML (20,000PPM)	FLUOROURACIL 50.0 MG/ML (50,000PPM)	VINCRIStINE SULFATE 1.0 MG/ML (1,000PPM)
DACARBAZINE 10.0 MG/ML (10,000PPM)	IFOSFAMIDE 50.0 MG/ML (50,000PPM)	FENTANYL CITRATE INJECTION 100mcg/2mL
DOXORUBICIN 2.0 MG/ML (2,000PPM)	MITOXANTRONE 2.0 MG/ML (2,000PPM)	