



# New Lab Start-Up PROGRAM



# We'll Set Up Your Lab from Start to Finish

## Program Overview

Government Scientific Source is here to make your dream laboratory become a reality. Our New Lab Start-Up Program provides exclusive support to setting up your new lab while potentially saving you significant start-up costs.

Our partnership with Laboratory Builders Inc, a leader in lab design and construction, allows us to assist in building a laboratory customized to your research needs. Our quick delivery lab program contains over 600 of the most commonly used lab products.

**Don't see something you need?** Not a problem, we can source it for you.

## Full-Service Offerings

GovSci's New Lab Start-Up program is more than just furniture in stock. Together with LBI, we can provide you with an experienced team of Furniture Specialists, Lab Designers, and Project Managers who will guide you through the entire process.

Our team offers:

- Clear and easy communication
- In-depth planning for labs of all disciplines
- Comprehensive project supervision from initial planning through final installation
- Design assistance, project management, and factory-trained installation
- An extensive inventory of top-quality furniture, equipment, devices and components, and space planning

# How it Works

Government Scientific's New Lab Start-Up Program gives you exclusive access to special offers that save you time and money. It streamlines your procurement process, making it much easier to secure everything you may need from us, your small-business one-stop shop.



## Step 1: Qualify

Our New Lab Start-Up Program is right for you if you are:

- Planning a new lab build
- Need retrofit or repairs to an existing lab
- Received funding for a new lab on your campus



## Step 2: Schedule

Schedule your consultation with a GovSci representative through our website or emailing [GSS-Sales@govsci.com](mailto:GSS-Sales@govsci.com). Once received, a rep will be in touch to get a better understanding of your unique project.

In addition, you can:

- Browse our quick delivery lab furnishings catalog of more than 600 lab products by category
- Shop from millions of lab products on govsci.com



## Step 3: Plan

- Work with your dedicated GovSci representative on building your lab from start to finish
- Identify items in the catalog
- Detail timelines and deliverables

## Schedule a Consultation ...



[govsci.com/promo/new-lab-start-up-program](http://govsci.com/promo/new-lab-start-up-program)



[GSS-Sales@govsci.com](mailto:GSS-Sales@govsci.com)



800.248.8030

# Your One-Stop Shop for Lab Supplies

We're here to support you even after your new laboratory build. Let Government Scientific Source be your reliable one-stop shop to advance your research.

GovSci's online marketplace includes millions of lab products from thousands of top-tier brands, including consumables, analytical instruments, chemicals, microscopy equipment, and a wide range of other scientific supplies. We simplify the procurement process with scientific expertise and dedicated customer service throughout your workflow, so that you have quick access to the products you need.



30+ years of government lab experience



75+ dedicated sales and customer service professionals



Millions of lab products from thousands of brands



Designated Small Business



GSA Advantage



Made in USA

## Contact Us ...



[govsci.com/contact](http://govsci.com/contact)




[GSS-Sales@govsci.com](mailto:GSS-Sales@govsci.com)



800.248.8030

# Categories



**NOTE:** The products listed in this brochure only feature a sample of the qualifying products available. Browse our full lab furnishings catalog to view the details of all products available for purchase: [govsci.laboutfit.com](http://govsci.laboutfit.com)

# Lab Casework ...

## Painted Steel Casework

### CornerStone *INSTOCK*® Standard Product Features:

- Inset style casework
- Modular design allows for simple field conversion from inset to overlay, door and drawer re-configurations and hardware component changes
- Durable construction provides SEFA 8 Metal compliant load-bearing capacity for years of dependable usage
- Drawer suspension is 100 lb. full extension
- *INSTOCK* “Pearl White” color is a urethane powder coating that provides superior chemical and abrasion resistance
- Products are SEFA 8 tested and approved by an independent third party



## Phenolic Casework

Phenolic resin casework is designed for rugged use in wet areas where corrosion or rust are of concern. Cabinets will not absorb odors, are easy to clean, and are virtually maintenance free. Any of our work surfaces, including stainless steel, can be combined with these cabinets to provide a complete laboratory solution.

Standard cabinet construction is from 1/2" thick chemically resistant phenolic panels in a full overlay style. All colors have a solid black core. Contact your Furniture specialist for assistance. Lead time is eight weeks.



## Stainless Steel Casework

### *INSTOCK*® Steel Laboratory Casework

The *INSTOCK* program is ready when you are and is designed to ship complete laboratories within 10 days.

Our *INSTOCK* program provides you an extensive selection of lab furnishings that are available in a matter of days.\*

All *INSTOCK* steel cabinets and chemical hoods have a powder coated “Pearl White” finish that is SEFA 8 compliant.\*\* A wide range of other colors and cabinet styles are available in six to eight weeks.



\*Most *INSTOCK*® items ship in 10 days. All items are subject to availability at the time of the order.

Large quantity orders may require longer lead times.

\*\*Contact us for a sample of Pearl White

Stainless steel furnishings are designed specifically for research and healthcare laboratories. Stainless steel counter tops and wash areas are also available for labs requiring sanitation and sterilization. Standard construction is a Type-304 stainless with a #4 finish. Our stainless steel cabinets have the same construction detail as our painted steel cabinets. Type-316 stainless cabinet construction is also available.

Contact your Furniture Specialist for more information. Lead time is six to eight weeks.



# Hoods ...

## Fume Hoods

### Hood Design Options to Fit Any Need with Unlimited Customization

ICI Scientific provides chemical fume hoods with exceptional safety, function, style and adaptability.



Our complete fume hood offering includes:

- True One Finger Sash Operation
- Easy Maintenance and Long Product Life Cycle
- Traditional Constant Volume
- ADA Accessible
- Full Access Walk-in
- Isotope
- CAV, VAV, and High Performance Hoods Available
- Distillation
- Demonstration
- Double-Sided
- Teaching
- Custom Enclosures
- Lifetime Warranty on Chain and Sprocket Sash Mechanism
- Perchloric Acid
- UL 1805 Listed
- Three-Year Warranty

## Ductless Fume Hoods

**Centrifugal Fan**  
Permanently lubricated fans are energy efficient with external rotor motor design that reduces operating costs. Extremely low noise and low vibration level.

**Sentinel Silver Microprocessor Control, Alarm System**  
Advanced microprocessor control supervises operation of all hood functions. Temperature-compensated air velocity sensor monitors airflow. 24-hour clock and blower run hour meter are standard.

**IN STOCK Quick Delivery Program**

**Powder Coat**  
Silver-Ion Impregnated powder coat. Inhibit microbial growth to improve safety.

**Robust Hood Construction**  
Key components including Fluorescent Lamps, Motor Capacitor, Electrical Harness, Electronix Ballast, and Switch Control are mounted outside the airstream and away from contaminated areas to permit easy service

**Sash Window**  
3 degree slope front allows easy access to work zone - enhancing ergonomics, eliminating operator fatigue and increasing productivity.

**Nanocarb Activated Carbon Filters**  
Optimized retention capacity. Diffusion technology to ensure even filter loading. Efficient perimeter clamping ensures no leakage as well as allows easy filter change procedure.

**Auto Purge Slots**  
Improves containment and operator protection by preventing accumulation of fumes in the work zone.

**Rear Walls**  
Electro galvanized steel wall for superior durability or tempered glass back wall.

**Work Top**  
1.5 mm, 16 gauge, type 304 stainless steel 4B finish. Easy to clean, corrosion and chemical resistant makes it an ideal worktop for experiments.

	Chemical Fume Containment	Filter Efficiency	Electrical Safety
Standards Compliance	ANSI/ASHRAE 110-1995, USA BS 7258, UK AFNOR NF X 15-203, France EN14175.3, Europe	BS 7989-2001, UK AFNOR NF X 15-211, France	UL-C-61010-1, USA CAN/CSA-22.2, No.61010-1 EN-61010-1, Europe IEC61010-1, Worldwide

# Cabinets ...

## Safety Cabinet

The Most Certified Biological Safety Cabinet in the World



Standards Compliance	Biosafety Cabinets	Air Quality	Filtration	Electrical Safety
NSF/ANSI 49	IEST-RP-CC001.3, USA IEST-RP-CC007, USA IEST-RP-CC034.1, USA	IEST-RP-CC001.3, USA IEST-RP-CC007, USA IEST-RP-CC034.1, USA	UL-C-61010A-1, USA IEC61010-1, Worldwide	

# Tables ...

## Patriot Tables

Our Patriot Table is a value driven, adjustable height, standard duty bench suitable for many basic lab applications. They are available with and without shelving uprights.



## Heavy Duty Steel Tables

Our Heavy Duty Steel Table is one of the most common tables used in laboratory settings and is typically located along the perimeter of a lab. It provides a sturdy and stable platform for your most critical work. It has a fully welded construction and is offered in fixed as well as adjustable heights.



# Shelving ...

## Stationary Wire and Mobile Wire Shelving

Shelving features:

- Chrome wire shelving is made of carbon-steel
- Stainless Steel Units are made of Type 304 Polished Stainless Steel
- Each unit includes: Four (4) posts, four (4) shelves, and split-sleeves to attach shelves to posts
- Open wire design minimizes dust accumulation and allows circulation for air
- Shelves can be adjusted at 1" intervals along entire height of the post
- 36"-54" long shelves have a maximum capacity of 800 lbs./each
- 60"-72" long shelves have a maximum capacity of 600 lbs./each
- Each mobile unit has a maximum capacity of 900 lbs., with uniformly distributed load
- Each stationary unit has a maximum capacity of 2,000 lbs., with uniformly distributed load
- Casters: Rubber casters with Poly housing, 5" swivel, two with brake



Stationary Wire Shelving



Mobile Wire Shelving

## Epoxy Coated Wire Shelving

Epoxy Coated Shelving Features:

- Each unit includes four (4) posts, four (4) shelves, and split-sleeves to attach shelves to posts
- Furnished in Metroseal 3 Epoxy
- Corrosion resistant finish for wet or high humidity environments
- Open wire design minimizes dust accumulation and allows free circulation for air
- Shelves can be adjusted at 1" intervals along entire height of the post
- 36"-54" long shelves have a maximum capacity of 800 lbs./each
- 60"-72" long shelves have a maximum capacity of 600 lbs./each
- Each mobile unit has a maximum capacity of 900 lbs., with uniformly distributed load
- Each stationary unit has a maximum capacity of 2,000 lbs., with uniformly distributed load
- Seven-year limited warranty against rust and corrosion
- Unit Casters: Rubber casters with Poly housing, 5" swivel, two with brake



Epoxy Coated Stationary Wire Shelving




Epoxy Coated Mobile Wire Shelving


# Workstations & Work Surfaces ...

## Envision Workstations

Our Envision Workstation provides plumbed and wired services that are in demand for today's self-contained & feature driven lab workstation. They can be placed in a pod with an overhead service panel or along the perimeter of a lab. An independent back-to-back design combined with an ambidextrous service layout make this workstation truly adaptable.



## Flexcore Laboratory Equipment Workstation



## Custom Manufactured Stainless Steel Work Surfaces

We offer custom laboratory work surfaces in stainless steel. These custom work surfaces are designed and engineered to suit your specific needs.

Sink and fixture locations will be sized and placed to suit your application. Lead time for stainless steel work surfaces is 7-8 weeks. Contact us for assistance.

30" Deep with 4" Rear Integral Splash	Unit
GSSSCUST	Per Linear Foot

All other configurations will be priced on a case by case basis.





# Chairs & Seating ...

## Customizable Seating: Laboratory and Industrial

Select from thousands of seating options, including chairs specifically designed for laboratory, clean room and ESD applications. Select the height, materials, features and colors that are most appropriate for your specific tasks or environment.

Ergonomics are critical, so our chairs and stools are designed to allow user to make adjustments to seats and backs and to add arms and foot rings where appropriate to provide optimal comfort and support.

Visit our chair configurator [gsschair.govsci.com](http://gsschair.govsci.com). Build the seat that is right for you.



Laboratory Seating



Industrial Seating

# Appliances & Accessories ...

## Faucets

### Mixing Faucets

- Mixing faucets are deck-mounted and have an option for a vacuum breaker
- They come with four-arm or blade handles
- They are chrome plated brass with an 8" swivel goose neck
- Valves are self-contained units with replaceable stainless steel seats



Part Number	Type
GSL412-8VB	Four-arm handle with vacuum breaker



Part Number	Type
GSL412-8	Four-arm handle without vacuum breaker



Part Number	Type
GSL412-8VB-BH	Blade handle with vacuum breaker



Part Number	Type
GSL412-8-BH	Blade handle without vacuum breaker

## Cold Water Faucets

- Cold water faucets are deck-mounted and have an option for a vacuum breaker.
- They come with a four arm or blade handle
- They are chrome plated brass with an 8" swivel goose neck
- Valves are self-contained units with replaceable stainless steel seats



Part Number	Type
GSL611-8VB	Four-arm handle with vacuum breaker

Part Number	Type
GSL611-8	Four-arm handle without vacuum breaker

Part Number	Type
GSL611-8VB-BH	Blade handle with vacuum breaker

Part Number	Type
GSL611-8-BH	Blade handle without vacuum breaker

## Drench Shower/Eyewash Stations

### Drench Shower/Eyewash

- Free-standing drench shower made of orange ABS plastic
- Ball valve operated by a rigid pull rod and easy grip handle
- Eyewash fixture with stainless steel bowl is provided with four gentle spray heads and a "Push to Operate" ball valve
- Eyewash and shower operated independently
- Overall height is 80"



Part Number	Type
GSSS909	Drench Shower/Eyewash

### Drench Shower

- Emergency shower head assembly is wall mounted
- It includes a 10" diameter orange ABS plastic shower head
- Operated by a rigid pull-rod and easy to grip handle



Part Number	Type
GSES643	Drench Shower

### ADA Drench Shower/Eyewash

- Free-standing drench shower made of orange ABS plastic
- Ball valve operated by a rigid pull rod and easy grip handle
- Eye wash fixture with stainless steel bowl is provided with four gentlespray heads and a "Push to Operate" ball valve
- Eyewash and shower operated independently
- Overall height is 80"
- Unit complies with ADA requirements for accessibility byhandicapped persons (maximum height and reach, minimum knee clearance and distance from obstructions)



Part Number	Type
GSSSBF909	Drench Shower

## Stainless Steel Sinks

### Drop-In Stainless Steel Sinks

- These stainless steel sinks are of a seamless die-drawn construction of 18 gauge, Type 304 stainless steel
- Interior surfaces are polished to a non-porous finish
- The underside is fully coated and insulated for sound deadening and reduced condensation
- Self-rimming top mount with stainless steel mounting channels
- The sink bottoms are pitched to a center drain
- Crumb strainer and tail piece included



Part Number	Dimensions (L x W x D)
GSS-1721-A-GR	21" x 17" x 7.5"
GSS-1925-A-GR	25" x 19" x 7.5"
GSSX-1925-A-GR	25" x 19" x 10.5"

### Undermount Stainless Steel

- These stainless steel sinks are of a seamless die-drawn construction of 18 gauge, Type 304 stainless steel
- Interior surfaces are polished to a non-porous finish
- The underside is fully coated and insulated for sound deadening and reduced condensation
- The sink bottoms are pitched to a center drain
- Crumb strainer and tail piece included
- Sink support required and ordered separately



Part Number	Dimensions (L x W x D)
GSUSX-1618-A	18" x 16" x 10.5"
GSUS-1824-A	24" x 18" x 7.5"
GSUSX-1824-A	24" x 18" x 10.5"
GSUSX-1830-A	30" x 18" x 10.5"

## Stainless Steel Garment Racks, Utility Cart, Security Units

### Garment Racks: Stainless Steel

- Electro-polished, Type 304 Stainless Steel hangers are included. Hangers are spaced 3" (7.6 cm) intervals
- Open wire construction promotes laminar air flow
- Unique hanger capture mechanism provides theft deterrence
- Easy to assemble
- Adjustable stainless steel triangular foot plates provided for floor mounting
- Upright unit is non-conductive



Upright



Wall Mount

Wall Mount	Upright	Dimensions W x D	No. of Hangers
GSGRW2436S	GSGRU2436S	36" x 24"	11
GSGRW2448S	GSGRU2448S	48" x 24"	15

### Gowning Bench: Stainless Steel

- Heavy-duty, solid top bench provides a sturdy and comfortable gowning fixture. All benches are 18" (45.7 cm) tall
- Top is manufactured from 14 gauge, Type 304 stainless steel
- Adjustable stainless steel foot plates provided for floor mounting
- Unit is non-conductive



Part Number	Dimensions W x D
GSGB1636S	36" x 16"
GSGB1648S	48" x 16"
GSGB1660S	60" x 16"
GSGB1672S	72" x 16"

### MetroMax Utility Cart

- Corrosion proof shelves, posts and Type 304 stainless steel handle
- 5" (12.7 cm) swivel casters, polyurethane wheels, poly housing and stainless steel axles
- Donut bumpers included with casters
- Furnished complete as a two shelf or three shelf unit



Part Number	Description	Dimensions W x D x H
GSMXUC1830-25	2 Shelves	18" x 30" x 39"
GSMXUC2436-25	2 Shelves	24" x 36" x 39"
GSMXUC1830-35	3 Shelves	18" x 30" x 39"
GSMXUC2436-35	3 Shelves	24" x 36" x 39"

### Security Units

- Durable chrome finish will not chip or flake
- Double doors each open 270° and can be secured along the sides of the unit
- Mobile units four (4) 5" polyurethane swivel casters with bumpers
- Unit comes with top and bottom shelves and 2 internal wire shelves
- Stationary Units are 67" tall. Mobile units are 68" tall



Stationary Unit

Part Number	Dimensions W x D x H
GSJSEC53	36" x 24" x 67"
GSJSEC55	48" x 24" x 67"
GSJSEC56	60" x 24" x 67"



Mobile Unit

Part Number	Dimensions W x D x H
GSJSEC53	36" x 24" x 67"
GSJSEC55	48" x 24" x 67"
GSJSEC56	60" x 24" x 67"

### Deep Ledge Polymer Utility Cart

- Specially designed with a 2 ¾" deep ledge to contain product and spills
- Polymer shelves are easy to clean and corrosion proof
- On units with three (3) shelves the center shelf adjusts at 1" increments
- Weight capacity is 400 lbs.
- Blue units have built-in Microban antimicrobial product protection



Grey	Black	Blue	Dimensions W x D x H
<b>2 Shelves</b>			
GSBC2030-2DG	GSBC2030-2DBL	GSBC2030-2DMB	21 ½" x 33 ¾" x 33 ¼"
GSBC2636-2DG	GSBC2636-2DBL	GSBC2636-2DMB	27" x 39 ½" x 33 ¼"
<b>3 Shelves</b>			
GSBC2030-3DG	GSBC2030-3DBL	GSBC2030-3DMB	21 ½" x 33 ¾" x 33 ¼"
GSBC2636-3DG	GSBC2636-3DBL	GSBC2636-3DMB	27" x 39 ½" x 33 ¼"

### Polymer Utility Cart

- Each cart comes with two labels for easy identification
- All-polymer shelves will never corrode
- Impact resistant shelves are designed to withstand daily abuse
- 7/16" (1.1 cm) shelf's edge retains small spills and prevents items from sliding off during transport
- Blue units have built-in Microban antimicrobial product protection



Grey	Black	Blue	Dimensions W x L x H
<b>2 Shelves</b>			
GSMY1627-24G	GSMY1627-24BL	GSMY1627-24BU	18 5/16" x 31 ½" x 35 ½"
GSMY2030-24G	GSMY2030-24BL	GSMY2030-24BU	23 7/16" x 34 3/8" x 35 ½"
GSMY2636-25G	GSMY2636-25BL	GSMY2636-25BU	27 11/16" x 40 ¼" x 36 7/8"
<b>3 Shelves</b>			
GSMY1627-34G	GSMY1627-34BL	GSMY1627-34BU	18 5/16" x 31 ½" x 35 ½"
GSMY1627-34G	GSMY2030-34BL	GSMY2030-34BU	23 7/16" x 34 3/8" x 35 ½"
GSMY2636-35G	GSMY2636-35BL	GSMY2636-35BU	27 11/16" x 40 ¼" x 36 7/8"

### Stainless Steel Work Tables with Bottom Shelf

- Seamless 14 gauge Type 304 stainless steel work surface and support structure
- Stationary work tables are 34" tall with stationary posts and leveling feet



Part Number	Dimensions D x L x H
GSWT305FS	30" x 48" x 34"
GSWT306FS	30" x 60" x 34"
GSWT307FS	30" x 72" x 34"
GSWT309FS	30" x 96" x 34"
GSWT366FS	36" x 60" x 34"
GSWT367FS	36" x 72" x 34"
GSWT369FS	36" x 96" x 34"
GSWT446FS	44" x 60" x 34"
GSWT447FS	44" x 72" x 34"
GSWT449FS	44" x 96" x 34"

### Stainless Steel Work Tables with Three-sided Frame

- Seamless 14 gauge Type 304 stainless steel work surface and support structure
- Stationary work tables are 34" tall with stationary posts and leveling feet



Part Number	Dimensions D x L x H
GSWT305US	30" x 48" x 34"
GSWT306US	30" x 60" x 34"
GSWT307US	30" x 72" x 34"
GSWT309US	30" x 96" x 34"
GSWT366US	36" x 60" x 34"
GSWT367US	36" x 72" x 34"
GSWT369US	36" x 96" x 34"
GSWT446US	44" x 60" x 34"
GSWT447US	44" x 72" x 34"
GSWT449US	44" x 96" x 34"

### Mobile Stainless Steel Work Tables

- Seamless 14 gauge Type 304 stainless steel work surface and support structure
- Mobile work tables are 34" tall with special posts and four (4) poly casters (two (2) swivel/two (2) brake)



Part Number	Dimensions D x L x H
GSMWT305US	30" x 48" x 34"
GSMWT306US	30" x 60" x 34"
GSMWT307US	30" x 72" x 34"
GSMWT309US	30" x 96" x 34"

### Balance Table

- Balance tables provide a stable and vibration resistant workspace for sensitive procedures
- Work top is 2" thick
- Top, sides and middle supports are all constructed from epoxy black resin material



Part Number	Dimensions W x D x H
GSBT3624	36" x 24" x 30"
GSBT4824	48" x 24" x 30"



**Distributed by Government Scientific Source.**

Contact Us: [GSS-Sales@govsci.com](mailto:GSS-Sales@govsci.com) • 800.248.8030 • [govsci.com](http://govsci.com)

© 2024 Government Scientific Source, Inc. All rights reserved.