

Enabling your next discovery



Save 25% on a dual-purpose CO₂-resistant shaker*

The Thermo Scientific™ Dual-Purpose CO₂-Resistant Orbital Shaker provides reliable around-the-clock operation ideally suited to keeping your cells alive and flourishing within your work environment. The CO₂-resistant shaker has the following features and benefits:

- Orbital shaking through magnetic power eliminates the need to replace worn belts
- Shaker has a small footprint, which is ideal for use in a CO₂ incubator
- Specially-treated mechanical components help protect exposure of up to 20% CO₂ and 95% humidity conditions
- External control box with LED display conveniently allows for adjustments without opening the chamber door

- Minimal heat-dissipation helps ensure conditions within the CO₂ chamber are not effected
- Stainless-steel construction protects against corrosive CO₂-gassing conditions
- An extensive selection of accessories fits most applications
- Dual-purpose—shaker can also be used outside of a CO₂ incubator for maximum value



Use promotional Cat. No. **88881103FLASH** when placing your order to get these savings.

Description	Promotional price	Cat. No.
CO ₂ Resistant Shaker Package	\$2,999.00	88881103FLASH
Package includes:		
CO ₂ Resistant Shaker		88881101
Universal platform		88881122
Clamp starter kit		88881133

Find out more at thermofisher.com/acceleratingscience

ThermoFisher
SCIENTIFIC

This product is intended for General Laboratory Use. It is the customer's responsibility to ensure that the performance of the product is suitable for customer's specific use or application. © 2021 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. **PR1721306 0521**

* All promotions valid for orders placed between July 1, 2021 to September 30, 2021 by residents of the United States and Canada. Packages and discounted add-on products must appear on same order to receive discount. "Savings" or "Discount" indicated in any promotion herein represent the maximum percentage difference between the aggregate price for individual purchase of each item contained in the promotional package at the current standard MSRP and the MSRP for the promotional price for each promotional package. Valid only in connection with purchase of promotional packages identified herein using the part numbers listed. All pricing information and values based on current list prices. Promotional products will be shipped with purchased products. Not valid in connection with any other offer or discount. Void where prohibited by law or company policy. See representative for additional information. Please note that the value of the specific offer(s) the customer may receive from the manufacturer under this program is a "discount or other reduction in price" to the customer under section 1128B(b)(3)(A) of the Social Security Act [42 U.S.C. 1320a-7b (b)(3) (a)]. Accordingly, the customer shall disclose this and any other discounts or other reduction in price received under this program under any state or federal program which provides cost or charge-based reimbursements to the customer for the products or services purchased under this program.