

Add a **splash of color** to the lab with a **SPECIAL OFFER up to 64% off** on the colorful array of Thermo Scientific™ Finnpiquette™ F2 Manual Pipettes*

Finnpiquette F2

Rugged construction and ease of maintenance for your more demanding applications



Rugged and durable

The Finnpiquette F2 is made of tough PVDF components that stand up to harsh chemicals and the damaging effects of UV light.

Fully-Autoclavable

Decontamination of the Finnpiquette F2 is simple because there's no need to disassemble/re-assemble the F2 for autoclaving, minimizing disruption and downtime.

Maintenance is Easy

Simply detach the tip cone for efficient daily maintenance or decontamination when not using an autoclave.



Thermo Scientific Finnpiquette F2 Manual Pipettes

Cat. No.	Description	Special Offer Price (USD)
Single channels		
4642010	Finnpiquette F2, 0.2-2 µl, µTip	\$180.00
4642020	Finnpiquette F2, 0.5-5 µl, µTip	\$180.00
4642030	Finnpiquette F2, 1-10 µl, µTip	\$180.00
4642060	Finnpiquette F2, 2-20 µl	\$180.00
4642130	Finnpiquette F2, 5-50 µl	\$180.00
4642070	Finnpiquette F2, 10-100 µl	\$180.00
4642080	Finnpiquette F2, 20-200 µl	\$180.00
4642090	Finnpiquette F2, 100-1000 µl	\$180.00
4642100	Finnpiquette F2, 0.5-5 ml	\$180.00
4642110	Finnpiquette F2, 1-10 ml	\$180.00
Multichannels		
4662000	Finnpiquette F2, 8-channel, 1-10 µl	\$330.00
4662010	Finnpiquette F2, 8-channel, 5-50 µl	\$330.00
4662020	Finnpiquette F2, 8-channel, 10-100 µl	\$330.00
4662030	Finnpiquette F2, 8-channel, 30-300 µl	\$330.00
4662040	Finnpiquette F2, 12-channel, 1-10 µl	\$400.00
4662050	Finnpiquette F2, 12-channel, 5-50 µl	\$400.00
4662060	Finnpiquette F2, 12-channel, 10-100 µl	\$400.00
4662070	Finnpiquette F2, 12-channel, 30-300 µl	\$400.00

For optimal results, accuracy, and precision, use Finntip Flex Tips with your Finnpiquette F2 Pipette

Multi-Channel 8	Multi-Channel 12	Single Channel	Pipette Description	Bulk	Rack Non Sterile	Rack Sterile	Filter Sterile
		4642010	Fp F2 0,2-2 µl	94060120	94060100	94060103	94056980
		4642020	Fp F2 0,5-5 µl	94060120	94060100	94060103	94056980
4662000	4662040	4642030	F2 1-10 µl	94060120	94060100	94060103	94056980
		4642060	F2 2-20 µl	94060320	94060310	94060313	94056510
4662010	4662050	4642130	F2 5-50 µl	94060320	94060310	94060313	94056520
4662020	4662060	4642070	F2 10-100 µl	94060320	94060310	94060313	94056520
		4642080	F2 20-200 µl	94060320	94060310	94060313	94056380
4662030	4662070	n/a	F2 30-300 µl	94060520	94060510	94060513	94056580
		4642090	F2 100-1000 µl	94060720	94060710	94060713	94056710
		4642100	F2 0,5-5 ml	9402030	9402070	9402073	94052550

*This special offer is open to customers in U.S. Valid for orders placed between July 1, 2020 and December 31, 2020. Please contact your sales representative for full details. Subject to dealer participation. Cannot be combined with other discounts or promotion. Offer void where prohibited, licensed or restricted by federal, state, provincial, or local laws or regulation or agency/institutional policy. Other restrictions may apply. U.S. customers acknowledge that this offer may include a discount or other price reduction that must be properly and accurately accounted for and reported by customer in accordance with all federal and state laws, including without limitation the federal anti-kickback law (42 C.F.R. § 1320a-7b(b)(3)(A)) and regulations thereunder (42 C.F.R. §1001.952(h)).

Find out more at thermofisher.com/finnpiquette

This product is intended for General Laboratory Use. It is the customer's responsibility to ensure that the performance of the product is suitable for customers' specific use or application © 2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. PFF2SPLASH 0620