

Enabling your
next discovery



Save 15% on select Nalgene beakers, tube racks and LDPE bottles*

Nalgene Narrow-Mouth LDPE Bottles with Closure

- Excellent chemical resistance to most acids, bases and alcohols make this a versatile bottle
- Can be used in freezer to -100°C (-148°F) for long-term storage
- Translucence allows easy viewing of liquid levels
- Linerless caps provide the ultimate in leakproof protection

Nalgene™ Unwire™ Test Tube Racks: Resmer™ Manufacturing Technology

- Smooth finish leaves fewer places for contaminants to lodge; no welds to break, no sharp edges.
- Endplates accept large labels for easy identification of samples
- Higher specific gravity than water so it stays submerged in a water bath and won't float
- Withstands most lab chemicals and temperatures up to 121°C for a variety of lab procedures
- Available in six bright colors to help keep samples organized
- No assembly is required; rack is ready to use

Nalgene Narrow-Mouth PPCO Bottles with Closure

- Linerless caps provide the ultimate in leakproof protection without the use of a liner that can wrinkle, cause leaks, or contaminate your reagents
- Autoclavable
- Made only from high-quality, laboratory-grade plastic materials for dependably low leachables and extractables

Nalgene PMP Low-Form Plastic Beakers

- Transparent PMP material allows for a clear view of beaker contents
- No-drip, single-spout design provides precision pouring control
- With ribs for easy stacking
- Silk screened with graduations, size code, resin code, catalog number, max. temperature and "No Flame" symbol for easy identification and safe use
- Autoclavable/Graduated/Transparent
- Meets ISO/DIN 7056 International Standards for Laboratory Plasticware Accuracy

Description	Promotional Cat. No.
Bottles	
Bottle, narrow-mouth, LDPE, 1 oz, 30 mL	2003-0001
Bottle, narrow-mouth, LDPE, 2 oz, 60 mL	2003-0002
Bottle, narrow-mouth, LDPE, 4 oz, 125 mL	2003-0004
Bottle, narrow-mouth, LDPE, 8 oz, 250 mL	2003-0008
Bottle, narrow-mouth, LDPE, 16 oz, 500 mL	2003-0016
Bottle, narrow-mouth, LDPE, 32 oz, 1,000 mL	2003-0032
Bottle, narrow-mouth, LDPE, 1/4 oz, 8 mL	2003-9025
Bottle, narrow-mouth, LDPE, 1/2 oz, 15 mL	2003-9050
Bottle, narrow-mouth, PP, 1 oz, 30 mL	2006-0001
Bottle, narrow-mouth, PP, 2 oz, 60 mL	2006-0002
Bottle, narrow-mouth, PP, 4 oz, 125 mL	2006-0004
Bottle, narrow-mouth, PP, 8 oz, 250 mL	2006-0008
Bottle, narrow-mouth, PP, 16 oz, 500 mL	2006-0016
Bottle, narrow-mouth, PP, 32 oz, 1000 mL	2006-0032
Bottle, narrow-mouth, PP, 1/4 oz, 8 mL	2006-9025
Bottle, narrow-mouth, PP, 1/2 oz, 15 mL	2006-9050
Bottle, narrow-mouth, PP, 1/8 oz, 4 mL	2006-9125
Narrow-mouth Economy Bottle, PP, 30 mL	2087-0001
Narrow-mouth Economy Bottle, PP, 60 mL	2087-0002
Narrow-mouth Economy Bottle, PP, 125 mL	2087-0004
Narrow-mouth Economy Bottle, PP, 250 mL	2087-0008
Bottle, wide-mouth, LDPE, 1 oz, 30 mL	2103-0001
Bottle, wide-mouth, LDPE, 2 oz, 60 mL	2103-0002
Bottle, wide-mouth, LDPE, 4 oz, 125 mL	2103-0004
Bottle, wide-mouth, LDPE, 8 oz, 250 mL	2103-0008
Bottle, wide-mouth, LDPE, 16 oz, 500 mL	2103-0016
Bottle, wide-mouth, LDPE, 32 oz, 1,000 mL	2103-0032
Bottle, wide-mouth, PP, 1 oz, 30 mL	2105-0001
Bottle, wide-mouth, PP, 2 oz, 60 mL	2105-0002
Bottle, wide-mouth, PP, 4 oz, 125 mL	2105-0004
Bottle, wide-mouth, PP, 8 oz, 250 mL	2105-0008
Bottle, wide-mouth, PP, 16 oz, 500 mL	2105-0016
Bottle, wide-mouth, PP, 32 oz, 1,000 mL	2105-0032
Bottle, wide-mouth, square, PP, 2 oz, 60 mL	2110-0002
Bottle, wide-mouth, square, PP, 6 oz, 175 mL	2110-0006
Bottle, wide-mouth, square, PP, 8 oz, 250 mL	2110-0008
Bottle, wide-mouth, square, PP, 16 oz, 500 mL	2110-0016
Bottle, wide-mouth, square, PP, 32 oz, 1,000 mL	2110-0032
Bottle, large, narrow-mouth, PP, 1/2 gal., 2 L	2203-0005
Bottle, large, narrow-mouth, PP, 1 gal., 4 L	2203-0010
Bottle, large, narrow-mouth, PP, 2 gal., 8 L	2203-0020

Description	Promotional Cat. No.
Bottle, amber, large, PP, 1/2 gal., 2 L	2204-0005
Bottle, amber, large, PP, 1 gal., 4 L	2204-0010
Bottle, amber, large, PP, 2 gal., 8 L	2204-0020
Bottle, large, narrow-mouth, PC, 500 mL	2205-0016
Bottle, large, narrow-mouth, PC, 1,000 mL	2205-0032
Bottle, large, narrow-mouth, PC, 4 L	DS2205-0010
Bottle, large, narrow-mouth, PC, 8 L	DS2205-0020
Bottle, large, narrow-mouth, PC, 2 L	DS2205-0210
Bottle, large, narrow-mouth, PC, 2.5 L	DS2205-0250
Beakers	
Nalgene PMP Low-Form Plastic Beakers, 4,000 mL	1203-4000
Nalgene PMP Low-Form Plastic Beakers, 50 mL	1203-0050
Nalgene PMP Low-Form Plastic Beakers, 100 mL	1203-0100
Nalgene PMP Low-Form Plastic Beakers, 150 mL	1203-0150
Nalgene PMP Low-Form Plastic Beakers, 250 mL	1203-0250
Nalgene PMP Low-Form Plastic Beakers, 400 mL	1203-0400
Nalgene PMP Low-Form Plastic Beakers, 600 mL	1203-0600
Nalgene PMP Low-Form Plastic Beakers, 1,000 mL	1203-1000
Nalgene PMP Low-Form Plastic Beakers, 2,000 mL	1203-2000
Unwire Test Tube Racks	
White 13 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0013
White 16 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0016
White 20 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0020
White 25 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0025
White 30 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0030
Orange 13 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0113
Orange 16 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0116
Orange 20 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0120
Orange 25 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0125
Orange 30 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0130
Yellow 13 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0213
Yellow 16 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0216
Yellow 20 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0220
Yellow 25 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0225
Yellow 30 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0230
Blue 13 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0313
Blue 16 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0316
Blue 20 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0320
Blue 25 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0325
Blue 30 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0330
Green 13 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0413
Green 16 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0416
Green 20 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0420
Green 25 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0425
Green 30 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0430
Red 13 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0513
Red 16 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0516
Red 20 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0520
Red 25 mm Nalgene Unwire Test Tube Rack, 1/8/pk/cs	5970-0525



Find out more at thermofisher.com/acceleratingscience

ThermoFisher
SCIENTIFIC

This product is intended for General Laboratory Use. It is the customer's responsibility to ensure that the performance of the product is suitable for customer's specific use or application. © 2021 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. **PR1850773 0521**

*All promotions valid for orders placed between July 1, 2021 to December 31, 2021 by residents of the United States (excluding Puerto Rico). Valid only in connection with purchase of promotional packages identified herein using the part numbers listed. All pricing information and values based on current list prices. Not valid in connection with any other offer or discount. Void where prohibited by law or company policy. See representative for additional information. Please note that the value of the specific offer(s) the customer may receive from the manufacturer under this program is a "discount or other reduction in price" to the customer under section 1128B(b)(3)(A) of the Social Security Act [42 U.S.C. 1320a-7b (b)(3)(a)]. Accordingly, the customer shall disclose this and any other discounts or other reduction in price received under this program under any state or federal program which provides cost or charge-based reimbursements to the customer for the products or services purchased under this program.