

Enabling your
next discovery



Save 15% on Nalgene tubing*



Nalgene Metric Non-phthalate PVC tubing

- Thermo Scientific™ Nalgene™ Metric Non-phthalate Clear Plastic PVC Tubing is manufactured from the highest-grade resin materials and performs like classic Nalgene 180 tubing but is free from orthophthalate DEHP plasticizers.
- Complies with FDA CFR 21 for food packaging, USP VI, Cytotoxicity, REACH, Latex-free, RoHS, and complies with California proposition 65.
- Resists hardening and discoloration from a wide range of laboratory chemicals and does not impart odor or taste making this tubing ideal for a wide range of non-medical, general lab, research and pharmaceutical applications.

Non-phthalate PVC vacuum tubing

- Thermo Scientific™ Nalgene™ Non-phthalate Clear Plastic PVC Vacuum Tubing is manufactured from the highest-grade resin materials and performs like classic Nalgene 180 tubing but is free from orthophthalate DEHP plasticizers.
- Complies with FDA CFR 21 for food packaging, USP VI, Cytotoxicity, REACH, Latex-free, RoHS, and complies with California proposition 65.
- Resists hardening and discoloration from a wide range of laboratory chemicals and does not impart odor or taste. Withstands full vacuum (29.9 in. Hg) at room temperature.

Ordering information

Product	Cat. No.
Nalgene™ Metric Non-phthalate PVC tubing	
PVC tubing, 25 m	8703-0406
PVC tubing, 25 m	8703-0407
PVC tubing, 25 m	8703-0508
PVC tubing, 25 m	8703-0609
PVC tubing, 25 m	8703-0610
PVC tubing, 25 m	8703-0710
PVC tubing, 25 m	8703-0811
PVC tubing, 25 m	8703-0812
PVC tubing, 25 m	8703-0913
PVC tubing, 25 m	8703-1014
PVC tubing, 25 m	8703-1216
PVC tubing, 25 m	8703-1520
Non-phthalate PVC vacuum tubing	
Nalgene Non-phthalate PVC vacuum tubing	8702-1210
Nalgene Non-phthalate PVC vacuum tubing	8702-0055
Nalgene Non-phthalate PVC vacuum tubing	8702-0065
Nalgene Non-phthalate PVC vacuum tubing	8702-1065
Nalgene Non-phthalate PVC vacuum tubing	8702-0145
Nalgene Non-phthalate PVC vacuum tubing	8702-1145
Nalgene Non-phthalate PVC vacuum tubing	8702-0210

Product	Cat. No.
Non-phthalate PVC tubing	
Nalgene Non-phthalate PVC tubing	8701-0010
Nalgene Non-phthalate PVC tubing	8701-9020
Nalgene Non-phthalate PVC tubing	8701-9030
Nalgene Non-phthalate PVC tubing	8701-4060
Nalgene Non-phthalate PVC tubing	8701-9090
Nalgene Non-phthalate PVC tubing	8701-4120
Nalgene Non-phthalate PVC tubing	8701-4140
Nalgene Non-phthalate PVC tubing	8701-0160
Nalgene Non-phthalate PVC tubing	8701-0180
Nalgene Non-phthalate PVC tubing	8701-0190
Nalgene Non-phthalate PVC tubing	8701-0200
Nalgene Non-phthalate PVC tubing	8701-0260
Nalgene Non-phthalate PVC tubing	8701-0310
Nalgene Non-phthalate PVC tubing	8701-0350
Nalgene Non-phthalate PVC tubing	8701-0020
Nalgene Non-phthalate PVC tubing	8701-0030
Nalgene Non-phthalate PVC tubing	8701-0060
Nalgene Non-phthalate PVC tubing	8701-0070
Nalgene Non-phthalate PVC tubing	8701-0080
Nalgene Non-phthalate PVC tubing	8701-0090
Nalgene Non-phthalate PVC tubing	8701-0100
Nalgene Non-phthalate PVC tubing	8701-0110
Nalgene Non-phthalate PVC tubing	8701-0120
Nalgene Non-phthalate PVC tubing	8701-0130
Nalgene Non-phthalate PVC tubing	8701-0140

Contact your sales representative:

Name:

Email:

Phone:

Find out more at thermofisher.com/acceleratingscience

ThermoFisher
SCIENTIFIC

This product is intended for General Laboratory Use. It is the customer's responsibility to ensure that the performance of the product is suitable for customer's specific use or application. © 2021 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. **PR1850725 0521**

*Offer open to customers in the US. Offer valid for orders received between July 1, 2021 and December 31, 2021, or until promotional supplies are depleted, whichever comes first. Please order using these special catalog numbers to take advantage of special pricing on products listed. Please ask your sales representative for details. Valid only in connection with the catalog numbers in this promotional brochure. "Discount", "savings" or similar terms represent difference between manufacturer's published list price and manufacturer's suggested price for promotional offer. Manufacturer reserves the right to discontinue or modify this promotion at any time. Customer acknowledges that this offer may include a discount or other price reduction that must be properly and accurately accounted for and reported by customer in accordance with all federal and state laws, including without limitation the federal antikickback law (42 C.F.R. § 1320a-7b(b)(3)(A)) and regulations thereunder (42 C.F.R. §1001.952(h)). Not valid in connection with any other offer or discount. Void where prohibited by law or company policy. Please note that items represented at no cost represent "in-kind discounts" and should be reported accordingly. Prices and specifications or products subject to change at any time.