



# HEI-PLATE



WHEN POWER MEETS SAFETY



[WWW.HEIDOLPHNA.COM](http://WWW.HEIDOLPHNA.COM)

## SAFE HEATING & STIRRING

The robust magnetic stirrers of the Hei-PLATE series were developed to optimally mix any viscosity of liquids from smooth to intense. They are perfectly suited for decomposing organic and inorganic substances. Our premium devices offer the highest operating safety and comfort standards. 600 W heating capacity, a strong magnetic field allows for easy stirring of high viscosity media, and two separate timers for rotation and heating. With a commitment to high quality products and dedicated customer service, your Heidolph purchase will be a worthwhile investment.



### KERA-DISK®

Our Kera-Disk® plate provides:

- Aluminum top plate allows for immediate heat transfer
- Ceramic coating makes the top plate chemical-resistant and scratch-proof.



### HEAT INDICATOR

Heating and hot surface light indicator.



### DIE-CAST HOUSING

The Hei-PLATE's fire-resistant aluminum die-cast housing protects against fumes, liquids, vapors and internal corrosion.



### TIMER

The time function allows you to define expiry times for the heating and rotation functions separately.

\*Only available on Hei-Connect



### SMART ROTATION

To protect you from splashes and burns, the rotation speed ramps up slowly.



### SAFETY CIRCUITS

To prevent laboratory fires, the Hei-PLATE includes two independent safety circuits which switch off heating in case of any overtemperature situation.

OFFER VALID THROUGH OCTOBER 2022

**Buy 3 get 1 free** - applicable with purchase of:

Hei-Tec- **036110550**

Hei-Tec w/ PT1000- **036110582**

---

### Heidolph North America

1235 N Mittel Blvd. Suite B, Wood Dale IL 601911  
224-265-9600 | [hello@heidolph.com](mailto:hello@heidolph.com) | [@heidolphna](https://www.heidolphna.com)