

Back to the Lab



Buy 3 cases of qualifying analytical filters, get 1 case at no additional cost*

Nalgene Sterile Analytical Filter Units and Filter Funnels

Recover and grow microorganisms for QC testing and analysis of water, food and beverage, raw materials and finished product with Thermo Scientific™ Nalgene™ Sterile Analytical Filter Units and Filter Funnels.

All devices feature an easily removable 47mm Cellulose Nitrate (CN) membrane and exceed ISO 7704 recovery testing.

- 0.2 µm is suitable for sterility testing
- 0.45 µm is certified for water quality testing and used for a wide range of beverage testing



Nalgene Analytical Filter Unit and Funnel

Membrane color	Gridded	Capacity (mL)	Pore size (µm)	Promo Cat. No.
Nalgene Analytical Test Filter Funnels, sterile				
White	No	100	0.2	145-0020
White	Yes	100	0.45	145-0045
White	No	250	0.2	145-2020
White	Yes	250	0.45	145-2045
Black/Grey	Yes	100	0.45	147-0045
Black/Grey	Yes	250	0.45	147-2045
Nalgene Analytical Test Vacuum Filter Units, sterile				
White	No	150	0.2	130-4020
White	Yes	150	0.45	130-4045

Maximum 3 redemptions per customer

*Valid in the United States (excluding Puerto Rico) and Canada on orders of 3 cases in a single purchase order. The case(s) at no additional charge will be of equal or lesser value than the purchased cases. Offer applies to orders received between April 9, 2020 and December 31, 2020 and only in connection with the products identified above. Maximum 3 redemptions per customer. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative or visit website for redemption details. Cannot be combined with other discounts or promotions. Offer void where prohibited, licensed or restricted by federal, state, provincial, or local laws or regulation or agency/institutional policy. Other restrictions may apply. Customer acknowledges that this offer may include a discount or other price reduction that must be properly and accurately accounted for and reported by customer in accordance with all federal and state laws, including without limitation the federal anti-kickback law (42 C.F.R. § 1320a-7(b)(3)(A)) and regulations thereunder (42 C.F.R. §1001.952(h)).

Find out more at thermofisher.com/backtothelab

This product is intended for General Laboratory Use. It is the customer's responsibility to ensure that the performance of the product is suitable for customer's specific use or application © 2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. PFNAFUBTTLICP 0820