

# SARTORIUS

## Special Entris<sup>®</sup> II Offer

### Best Value in its Class for Basic Weighing Tasks

Every Sartorius balance offers quality, value, and consistency

No matter what you're weighing, the new Entris<sup>®</sup> II balance is always the right choice. Offering unrivalled value and backed by almost 150 years of German engineering expertise, the Entris<sup>®</sup> II comes in two product lines, so you can find the balance that meets your specific weighing needs.

Our Entris<sup>®</sup> II Essential line is the only balance in its class offering isoCAL, LED touch technology, and 12 built-in applications at a budget price. With over 40 models from which to choose, you will surely find a model to meet your basic weighing needs.



Get a Midi Plus<sup>™</sup> Controller  
at no additional cost!

# Special Midi Plus™ Controller Offer

Buy select models of Entris® II BCE and get a Midi Plus™ Controller at no cost.

Midi Plus™ is a lightweight, cordless pipetting controller with excellent ergonomics, enabling you to aspirate from bottles and tubes without holding your arm and hand in an uncomfortable raised position.

- Precise finger tip control of variable speed ranges
- Fits all commonly used 1 - 100 ml glass or plastic pipettes
- Hydrophobic autoclavable filter prevents from overaspiration
- Low battery warning



Special offer valid only for the US from April 1, 2021 through September 30, 2021. Redemptions must be submitted by October 15, 2021. This offer cannot be used in combination with other special offers. One no charge Midi Plus per customer location.

Model Number    Weighing capacity (g)    Readability (mg)    List Price

#### Entris II Basic Essential Analytical Balance with Internal Calibration

BCE124I-1S	120	0.1	\$2,600.00
BCE224I-1S	220	0.1	\$2,910.00
BCE64I-1S	60	0.1	\$2,290.00

#### Entris II Basic Essential Precision Balance with Internal Calibration

BCE1202I-1S	1200	10	\$1,560.00
BCE2201I-1S	2200	100	\$1,455.00
BCE2202I-1S	2200	10	\$1,665.00
BCE223I-1S	220	1	\$1,455.00
BCE3202I-1S	3200	10	\$1,770.00
BCE323I-1S	320	1	\$1,665.00
BCE4202I-1S	4200	10	\$1,975.00
BCE423I-1S	420	1	\$2,080.00
BCE5201I-1S	5200	100	\$1,560.00
BCE6202I-1S	6200	10	\$2,495.00
BCE6200I-1S	6200	1000	\$1,455.00
BCE622I-1S	620	10	\$1,350.00
BCE623I-1S	620	1	\$2,290.00
BCE653I-1S	650	1	\$2,495.00
BCE8201I-1S	8200	100	\$1,770.00
BCE8200I-1S	8200	1000	\$1,560.00
BCE822I-1S	820	10	\$1,455.00

#### Get a Midi Plus at no cost

710931

## Featured & Benefits

Easy clean for operational efficiency and durability

- High chemical resistance
  - ensured using parts made from hard-wearing polybutylene terephthalate (PBT), stainless steel and glass
- Prevents cross contamination
  - with wipe clean design and easy-to-remove parts

Effective draft shield

- Minimizes weighing errors caused by electrostatically charged samples using specially coated glass parts

isoCAL (internal calibration and adjustment) for total assurance of accurate weighing results

- Optimal accuracy and operating convenience
  - using fully automated internal temperature and time-controlled calibration and adjustment feature, unique in this weighing class
- Assures (SOP) compliant operation
  - with self-notification if calibration is outside the normal range



Please contact your Sales Representative or any of our authorized Partners for more information

Find out more:  
[sartorius.com](http://sartorius.com)