

Enabling your next discovery



Save when you buy Thermo Scientific™ Precision™ General Purpose Water Baths and thermal beads*

Thermal beads offer a great alternative to water

Thermo Scientific™ Thermal Beads are small, nonuniform metal beads composed of dry, metallic thermal alloy, which are designed to replace water in a water bath. They conveniently heat samples with noncritical temperature requirements because of their size, shape, and density, which enable vessels of any size or shape that fit into the bath to be held in place without floating. Using thermal beads also helps reduce the risk of contamination because there is no liquid to circulate contaminants nor to promote biogrowth.

- Easy use and maintenance
- Helps reduce risk of contamination
- Designed to replace water in a water bath



Product bundle	Promotion price	Cat. No.
Precision General Purpose Water Bath, 2 L, and 4 L of Thermal Beads	\$1,502.97	TSGP02BDL
Precision General Purpose Water Bath, 5 L, and 4 L of Thermal Beads	\$1,646.36	TSGP05BDL
Precision General Purpose Water Bath, 10 L, and 8 L of Thermal Beads	\$2,796.66	TSGP10BDL
Precision General Purpose Water Bath, 20 L, and 16 L of Thermal Beads	\$4,284.14	TSGP20BDL
Precision General Purpose Water Bath, 28 L, and 24 L of Thermal Beads	\$6,593.53	TSGP28BDL
Precision General Purpose Water Bath, 5 L and 10 L dual, and 12 L of Thermal Beads	\$3,848.76	TSGP15DBDL

Enter a single Cat. No. to get the bundle of contents listed on each line.

Find out more at thermofisher.com/acceleratingscience

ThermoFisher
SCIENTIFIC

This product is intended for General Laboratory Use. It is the customer's responsibility to ensure that the performance of the product is suitable for customer's specific use or application. © 2021 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. **PR1739658 0521**

* Offers apply to new qualifying purchases between June 1, 2021, and December 31, 2021 in the United States and Canada. For more details regarding this promotion contact your local sales representative. Cannot be combined with other discounts or promotions. Offer void where prohibited, licensed or restricted by federal, state, provincial, or local laws or regulation or agency/institutional policy. Other restrictions may apply. Customer acknowledges that these offers may include a discount or other price reduction that must be properly and accurately accounted for and reported by customer in accordance with all federal and state laws, including without limitation the federal anti-kickback law (42 C.F.R. §1320a-7b(3)(A)) and regulations thereunder (42 C.F.R. §1001.952(h)).